Doro 9350

Pill dispenser

ENGLISH

In this device we have integrated a radio transceiver to communicate wirelessly between device and carephone. In order for the device to send alarms through the carephone, you need to pair them.

How to pair the radio transceiver with the carephone

1. **On the pill dispenser;** press and hold all three buttons until TEST is shown on the display, then press button 3 to initiate full alarm sequence.
2. **On the carephone;** press and hold the green button for four seconds until a beep sounds to start pairing mode. Release the button.
3. **On the carephone;** press the **alarm button** to assign the new device. A new confirmation beep is heard.
4. **On the pill dispenser;** press and hold button 1 to leave TEST mode.

For more information on pairing the carephone with the radio transceiver, refer to the manual of the carephone.

How to test the alarm after pairing

On the pill dispenser; press and hold all three buttons until TEST is shown on the display, then press button 3 to initiate full alarm sequence. Press and hold button 1 to leave TEST mode. Check that the alarm is received at the Alarm Receiving Center.

This is how it works

1. When an alarm is triggered on the device, the radio transceiver sends the alarm to the carephone.
2. The carephone passes the alarm to the alarm receiving centre or another receiver.

Battery warning

The device automatically alarms at low battery. Replace the battery as soon as possible after low battery level has been indicated.

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©2017 Doro AB. All rights reserved.